



Urgent Field Safety Notice
N. 02/2019

Product Name	List Number (LN)	Lot Number	Expiration Date
Lithium	8L25-30	n/a	n/a

Date: June 28th, 2019

Details on affected devices:

The purpose of this communication is to inform of an additional SmartWash parameter for the MULTIGENT Lithium reagent, List Number (LN) 8L25-30 to prevent the potential for carryover from the ARCHITECT Lactate Dehydrogenase (LDH) reagent (LN 2P56). Please review the information carefully and follow the mandatory actions.

Description of the problem:

Carryover may be observed between the MULTIGENT Lithium reagent and the ARCHITECT LDH reagent on improperly maintained systems due to the presence of lithium lactate in LDH reagent. As a result, falsely elevated lithium patient results may be generated. To minimize the potential for carryover, an additional SmartWash has been implemented for the MULTIGENT Lithium assay.

Patient Impact:

There is a potential to generate falsely elevated lithium patient results. Therefore, the following actions are mandatory.

Actions to be taken:

- If your laboratory **does not** have the LDH reagent installed on your system, no action is required.
- If the LDH reagent is installed on your instrument, please follow the instructions below:
 1. Update the Lithium (LN 8L25-30) assay parameters to add the new Smart Wash. For additional information, see the ARCHITECT System Operations Manual, Section 2, *Configure the SmartWash settings (c System)*.

Configure assay parameters – SmartWash				
O General	O Calibration	• SmanWash	O Results	O Interretation
Ass : Lith				
GOMponent	REAGENT/ASSAY	WASH	Volume	Replicates
R1	LDH00	Water		!
SampleProbe•		Water		
•sample Probe <i>Sample wash protocol</i> is Maximumwash .				

2. For information on maintaining optimal system performance, refer to your ARCHITECT System Operations Manual.
- Please review the content of this communication with your Medical Director and retain this letter for your laboratory records.



Transmission of this Field Safety Notice:

This notice needs to be passed on all those who need to be aware within your organization or to any organization/individuals where the potentially affected devices have been transferred.

Contact reference person:

If you or any of the health care providers you serve have any questions regarding this information, please contact your local area Customer Service.

Best regards,

...