

Rev 1: September 2018
FSN Ref: FSN-2021-0002
Date: 20 May 2021

FSCA Ref: FSN-2021-0002

Urgent Field Safety Notice
RapID NF Plus System, R8311005, Batch 158548

For Attention of*:

Contact details of local representative (name, e-mail, telephone, address etc.)*

E.mail : mbd.vigilance@thermofisher.com

Telephone: +44(0) 1256 841144

Fax: +44(0) 1256 479525

RapID NF Plus System, R8311005, Batch 158548

1. Information on Affected Devices*	
1.	1. Device Type(s)* IVD
1.	2. Commercial name(s) RapID NF Plus System
1.	3. Unique Device Identifier(s) (UDI-DI) 00848838058158
1.	4. Primary clinical purpose of device(s)* Remel RapID™ NF Plus System is a qualitative micromethod employing conventional and chromogenic substrates for the identification of medically important glucose non-fermenting, Gram-negative bacteria and other select glucose-fermenting, Gram-negative bacteria not belonging to the family Enterobacteriaceae, which have been isolated from human clinical specimens. A complete listing of the organisms addressed by the RapID NF Plus System is provided in the RapID NF Plus Differential Chart (found in the IFU).
1.	5. Device Model/Catalogue/part number(s)* R8311005
1.	6. Software version N/A
1.	7. Affected serial or lot number range 158548
1.	8. Associated devices N/A


2. Reason for Field Safety Corrective Action (FSCA)*	
2.	1. Description of the product problem* ATCC 19606 (<i>Acinetobacter baumannii</i> ATCC® 19606), ATCC 13253 (<i>Elizabethkingia meningoseptica</i> ATCC® 13253) and blank (NF reagent) gave a positive reaction where it should have given a negative reaction within the NO ₃ well of the panel.
2.	2. Hazard giving rise to the FSCA* There should be no immediate or long-term health consequences from using this product. The determination of nitrate in the affected species are not the sole determinant for identification of these species. There are some strains of both <i>A. baumannii</i> and <i>E meningosepticum</i> that are positive for NO ₃ , so the entire range of biochemical tests should be considered in the identification of clinical specimens. In this context of a single false positive test, the clinical risk should be considered negligible.
2.	3. Probability of problem arising High

Rev 1: September 2018

FSN Ref: FSN-2021-0002

FSCA Ref: FSN-2021-0002

2.	<p data-bbox="325 241 767 275">4. Predicted risk to patient/users</p> <p data-bbox="277 280 938 313">The clinical risk to patients is considered as minor.</p> <p data-bbox="277 344 1382 448">The panels are used for confirmation of the organism and when the panel results are entered into the supportive ERIC database the presumptive result is confirmed, with minimal impact of the overall result despite the NO₃ well failure.</p> <p data-bbox="277 450 1305 553">The probability of the identification when using ATCC 13253 remains as <i>Elizabethkingia meningoseptica</i> regardless if a positive or negative NO₃ result is inputted.</p> <p data-bbox="277 555 1382 714">However, for ATCC 19606 for the identification of <i>Acinetobacter baumannii</i> a higher probability of identification as <i>Burkholderia cepacia</i> is received. Here the ERIC database states that the frequency value for the first choice is not within acceptable limits for identification, and instructs the user to recheck coding, test interpretation, Gram stain morphology and test procedures, re-isolation, and reinput results.</p> <p data-bbox="277 716 1382 882">With this <i>Acinetobacter</i> strain discrepancy, the isolate is first reported as a non-fermenting Gram negative bacilli and if clinically important in the specimen then antimicrobial susceptibilities are performed. Patient will already be placed on broad spectrum agents. The clinical risk of the initial "error" in identification is minor with no delay in treatment.</p> <p data-bbox="277 884 1142 918">The patient treatment will remain unchanged regardless of results.</p>
2.	<p data-bbox="325 949 1054 983">5. Further information to help characterise the problem</p> <p data-bbox="373 985 427 1016">N/A</p>
2.	<p data-bbox="325 1016 655 1050">6. Background on Issue</p> <p data-bbox="277 1055 1347 1122">One customer complaint had been received. The complaint was confirmed as the transferred retained sample replicated the issue.</p> <p data-bbox="277 1124 1409 1191">This product was manufactured at a Thermo Fisher Scientific manufacturing site which is no longer in existence.</p>
2.	<p data-bbox="325 1225 836 1258">7. Other information relevant to FSCA</p> <p data-bbox="277 1261 1382 1294">Lot. 158548 was manufactured in 11-March-2020 with the expiry of 03-August-2021.</p>

4. General Information*		
4.	1. FSN Type*	New
4.	2. For updated FSN, reference number and date of previous FSN	N/A
4.	3. For Updated FSN, key new information as follows:	
	N/A	
4.	4. Further advice or information already expected in follow-up FSN? *	Not planned yet
4	5. If follow-up FSN expected, what is the further advice expected to relate to:	
	N/A	
4	6. Anticipated timescale for follow-up FSN	N/A
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Thermo Fisher Scientific
	b. Address	Clipper Boulevard West, Cross ways industrial estate, Dartford, Kent. DA2 6PT
	c. Website address	www.thermofisher.com
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
4.	9. List of attachments/appendices:	Customer response form
4.	10. Name Vice President, Quality and Regulatory, MBD
	Signature	

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.*</p>

Rev 1: September 2018
FSN Ref: FSN-2021-0002

FSCA Ref: FSN-2021-0002

Customer Reply Form

1. Field Safety Notice (FSN) information			
FSN Reference number*	FSN-2021-0002		
FSN Date*	20 May 2021		
Product/ Device name*	Remel RapID™ NF Plus System		
Product Code(s)	R8311005		
Batch/Serial Number (s)	158548		
2. Customer Details			
Account Number			
Organisation Name*			
Organisation Address*			
Department/Unit			
Shipping address if different to above			
Contact Name*			
Title or Function			
Telephone number*			
Email*			
3. Customer action undertaken on behalf of Healthcare Organisation			
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.		
<input type="checkbox"/>	I performed all actions requested by the FSN.		
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.		
<input type="checkbox"/>	I have returned affected devices - enter number of devices returned and date complete or N/A	Qty:	Lot/Serial Number: Date Returned (DD/MM/YY)
		Comments:	
<input type="checkbox"/>	I have destroyed affected devices – enter number destroyed and date complete. (EDIT WHEN NECESSARY) or N/A	Qty:	Lot/Serial Number: Date Returned (DD/MM/YY)
		Qty	Credit <input type="checkbox"/> Replacement <input type="checkbox"/>
		Comments:	
<input type="checkbox"/>	No affected devices are available for return/ destruction		
<input type="checkbox"/>	Other Action (Define):		
<input type="checkbox"/>	I do not have any affected devices.		
<input type="checkbox"/>	I have a query please contact me (e.g. need for replacement of the product).		
Print Name*			
Signature*			
Date*			
4. Return acknowledgement to sender			
Email	MBD.vigilance@thermofisher.com		
Telephone Number & Fax	Tel : +44(0) 1256 841144 Fax :+44(0) 1256 479525		
Postal Address			
Deadline for returning the reply form*	24 June 2021		

Mandatory fields are marked with *

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN. Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.