

URGENT FIELD SAFETY NOTICE
fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution
Multiple Issues related to Software Anomalies
for Field Safety Corrective Action FSCA-21-002

30 November 2021

FSN Ref: FSCA-21-002_FSCA-21-003-FSN-2

Attention: Distributors and users of the fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution ventilators

Dear Customer,

The purpose of this communication is to inform you that a software upgrade (software version 5.2.1) and revised Instructions for Use associated with FSCA-21-002, as communicated in the original Field Safety Notice (FSN) FSCA-21-002_FSCA-21-003-FSN-1, are now available for the fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution ventilators.

Details of affected devices

Affected versions of fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution
 (*Issue # as referenced in the FSN FSCA-21-002_FSCA-21-003-FSN-1)

Issue #*	Issue	Affected devices		
		fabian HFO	fabian +nCPAP evolution	fabian Therapy evolution
1	Interruption of High Frequency Oscillation (HFO) in HFO Ventilation mode	112001 113001	Not affected	Not affected
2	Presence of incorrect display of external bias flow	113001	Not affected	Not affected
3	Presence of no alarm on endotracheal tube (ETT) disconnection	112001 113001	Not affected	Not affected
4	Global Alarms Off function becomes enabled during ventilation	112001, 113001	122001	121001
5	Graphical User Interface (GUI) freeze	111001 111001.01 112001 113001	122001	121001
6	Pressure delivery is below specification with Infant Flow LP circuits - fabian Therapy evolution, fabian +nCPAP evolution, and fabian HFO	111001 111001.01 112001 113001	122001	121001

NOTE: Corrections to the Field Safety Notice FSCA-21-002_FSCA-21-003-FSN-1 regarding affected devices:

- The scope (devices affected) for Issue 5 was incomplete in FSN FSCA-21-002_FSCA-21-003-FSN-1. Issue 5 affects all fabian devices (fabian HFO, +nCPAP evolution, and Therapy evolution). Refer to the following table for the corrections.

Page #	Incorrect text	Correct Text
3	Issue 5: Graphical User Interface (GUI) freeze / Application error and potential shut-down – fabian HFO	Issue 5: Graphical User Interface (GUI) freeze / Application error and potential shut-down – fabian HFO, fabian +nCPAP, and fabian Therapy evolution
7	Issue 5: Graphical User Interface (GUI) freeze when using etCO2 module – fabian HFO with etCO2 module	Issue 5: Graphical User Interface (GUI) freeze / Application error and potential shut-down – fabian HFO, fabian +nCPAP, and fabian Therapy evolution

- The devices with reference numbers 121012 and 122012 that were listed on FSN FSCA-21-002_FSCA-21-003-FSN-1 are not in the scope of this FSCA and will be addressed separately.

Once the new software update (**fabian Software Release Package 5.2.1**) is installed, you should use the device according to the updated Instructions for Use provided by the distributor or service partner.

Note: fabian devices that have not yet had FSCA-18-004 or FSCA-20-001 implemented can be updated directly to the new software version 5.2.1. Distributors should refer to the technical bulletin *Technical Bulletin TB-0035 Release of Software Version 5.2.1* for further information about the software update strategy.

Next Steps

The Acutronic Distribution partner / authorized technical service engineer will inform end users about the new software via an *End User Release Note* and make the necessary arrangements to install the software on the affected device(s).

Actions to be taken by distributors / authorized technical service partners

- Log into your client account on the Acutronic client website as described in the *Technical Bulletin TB-0035 Release of Software Version 5.2.1*. (If you do not have an account, you need to register.)
- Download the Instructions for Use for the affected device(s) from the client FSCA download area.
- Download the relevant Software Release Package 5.2.1 from the Acutronic client server as described in the technical bulletin *Technical Bulletin TB-0035 Release of Software Version 5.2.1*. Choices are **fabian HFO**, **fabian Therapy evolution**, and **fabian +nCPAP evolution**.
- Check the contents of the download. Software Release Package 5.2.1 contains the following:
 - Release note
 - Technical Release Note
 - End User Release Note
 - PIC package for programmers
 - USB package
 - Software update description
 - Test instructions
 - fabian Field Safety Corrective Action - FSCA-21-002 Completion Data & Verification Record form

- Inform immediately the end users of the fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution ventilators in scope of this FSCA about the fabian Software Release 5.2.1 by sending them this *FSN*, the *End User Release Note* and the relevant update of the *Instructions for Use* for Software Release 5.2.1.
- Install the software upgrade according to the upgrade instructions.
- Perform calibration and testing according to the test instructions.
- Fill out a fabian Field Safety Corrective Action - FSCA-21-002 Completion Data & Verification Record form for each device successfully upgraded to version 5.2.1, and return it using the following email address:

GMB-AMS-FSCAresponsecentre@vyaire.com

Actions to be taken by the end users

- Make sure that across the healthcare facility, this *FSN*, the *End User Release Note* and the *Instructions for Use* for Software Release 5.2.1 are made available immediately to any potential user of the fabian HFO, fabian +nCPAP evolution, and fabian Therapy evolution ventilators in the scope of this FSCA.
- Make sure that all potential users are adequately trained according to local training protocols.
- If you have any questions regarding installation of the software, please refer to your Acutronic/Vyair Distribution/authorized technical service partner or Acutronic/Vyair Sales Representative, as appropriate.

Contact information

For questions, concerns or any event that reasonably suggests is related to the subject of this FSCA, please email GMB-AMS-FSCAresponsecentre@vyaire.com.

The undersigned confirms that this notice has been notified to the appropriate Regulatory Agencies.

Sincerely,

...