Urgente veiligheidskennisgeving

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A large grey rectangular box redacting the name and contact information of the Medtronic representative.

Medtronic
Engineering the extraordinary

Bijlagen:

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Appendix A - IFU summarized anticipated updates for OUS (text may change)

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The [Controller Fault] alarm indicates that a controller malfunction may have occurred, but the pump is still working. The [Controller Fault] alarm is designed to also be a signal that the controller has reached the end of its expected useful life.

Patients should bring all controllers, power sources (batteries, AC and DC adapters), and alarm adapters to clinic visits. During clinic visits, the healthcare provider or clinician should ensure that the controller is within the expected useful life:

- If the primary and back-up controllers are within expected useful life (2 years from when they were provided to the patient), then continue to routinely inspect the controller connections and their pins.
- If the back-up controller has reached the end of its expected useful life (2 years from when it was provided to the patient), take it out of service and replace it with a new controller.
- If the primary controller has reached the end of its expected useful life (2 years from when it was provided to the patient), download the log files and send them to HeartWare for analysis:
 - If the log file analysis indicates that there are [Controller Fault] alarms after the end of the expected useful life, refer to the [Controller Fault] alarm instructions.
 - If the log file analysis indicates that there have been no [Controller Fault] alarms after the end of the expected useful life, the clinician should decide whether the patient can continue to use the controller:
 - In the future, there may be a [Controller Fault] alarm if the internal battery reaches its end of life. This [Controller Fault] alarm will be an active alarm that does not resolve. The risk associated with internal battery end of life is that the controller may not sound the [No Power] alarm when both power sources are disconnected. The clinician should assess on an individual basis if the risk associated with the internal battery end of life outweighs the risk associated with performing a controller exchange. If a controller exchange is deemed necessary to address a [Controller Fault] alarm due to internal battery end of life, program a new controller instead of using the back-up controller.

Factors that should be considered before performing a controller exchange include, but are not limited to:

- Consider whether the patient can tolerate a pump stop during the controller exchange.
- Patient and caregiver understanding/compliance of alarm responses and power source management.
- Distance/time it will take for the patient to reach the hospital/clinic for support.
- Length of time the patient is expected to remain on therapy.

If a controller exchange is deemed necessary, confirm that the remaining useful life of the back-up controller exceeds the length of time the patient is expected to remain on therapy. The expected useful life of the primary and back-up controller is 2 years from when it was provided to the patient.

A [Controller Fault] alarm is designed to occur after the end of the controller's expected useful life if there is an internal battery end of life. This will be indicated in the log files analysis. This [Controller Fault] alarm will be an

Medtronic

active alarm that does not resolve. The alarm will continue to sound an intermittent beep unless it is permanently silenced by the clinician using a monitor. The risk associated with internal battery end of life is that the controller may not sound the [No Power] alarm when both power sources are disconnected. Thus, the clinician should assess on an individual basis if the risk associated with the internal battery end of life outweighs the risk associated with performing a controller exchange. If a controller exchange is deemed necessary to address a [Controller Fault] alarm due to internal battery end of life, program a new controller instead of using the back-up controller.

WARNING: A [Controller Fault] alarm may occur after the end of the controller's expected useful life. If the alarm is due to the controller's internal battery end of life, the controller may not sound the [No Power] alarm when both power sources are disconnected.

WARNING: Always assess if the patient can tolerate a pump stop before exchanging the controller. The pump will stop for the duration of a controller exchange.

The [Controller Fault] alarm audio can be permanently muted but this cannot be done by the patient. It requires a monitor. The patient is only able to mute the alarm for 5 min or 1 hr. Muting the alarm does not resolve the [Controller Fault] alarm.

Recommendations for back-up controller:

Once a week: Instruct the patient to inspect all 4 connections and their pins for dirt or debris. The patient should not attempt to clean the controller connectors. The patient should be instructed to contact their clinician if they notice the connectors contain dirt or debris. Exterior surfaces of the controller should be cleaned using a clean cloth. A damp cloth may be used but a wet cloth should not.

Recommendations for the HVAD batteries:

Ensure that the batteries are within their expected useful life:

- Download the controller log files to determine the number of times the battery has been charged and discharged. The expected useful life is 500 charge and discharge cycles; if the use of the batteries is rotated, this should provide patient support for 1 year. Batteries that reach the end of their expected useful life should be taken out of service and replaced. If a battery lasts less than 2 hours after being fully charged, it should be taken out of service and replaced.

Recommendations for all HVAD system components:

Inspect all HVAD system components for signs of damage. Damaged equipment should be reported to Medtronic and replaced.

The HVAD peripherals and accessories were designed and tested to function for the following periods starting from when it was provided to the patient:

HVAD system peripherals or accessories	Expected useful life (starting from when it was provided to the patient)
--	---

Medtronic

HVAD controller (primary and back-up)	2 years
HeartWare batteries	500 charge and discharge cycles
Controller AC adapter	6 months
Controller DC adapter	1 year
Alarm adapter	1 year
Monitor data cable	1 year
HeartWare battery charger	1 year
Carry cases (HeartWare shoulder pack, waist pack, convertible patient pack)	1 year

Medtronic

Appendix B - Patient Manual (PM) summarized updates

A consolidated summary of the language being added to the US PM is provided below. Note that this information may be provided within various sections of the PM as appropriate

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The [Controller Fault] alarm indicates that a controller malfunction may have occurred, but the pump is still working. The [Controller Fault] alarm is designed to also be a signal that the controller has reached the end of its expected useful life.

WARNING: A [Controller Fault] alarm may occur after the end of the controller's expected useful life. Do not exchange the controller. Always call your clinician for appropriate action. The controller may not sound the [No Power] alarm when both power sources are disconnected.

WARNING: If a controller has reached the end of its expected useful life, do not exchange your controller. Always call your clinician for appropriate action.

Recommendations for back-up controller:

Once a week, inspect all 4 connections and their pins on the controller for dirt or debris. Do not attempt to clean the controller connectors. If any dirt or debris is found, report the condition to your clinician.

Recommendations for the HVAD batteries:

Ensure that the batteries are within their expected useful life:

- Each fully charged battery provides approximately 4 to 7 hours of use for normal activities such as reading or watching TV. The battery may last for less time as your activity level increases. Similar to the battery in a cell phone (or mobile phone), the HeartWare batteries lose charge over time. If a battery lasts less than 2 hours after being fully charged, take it out of service and contact your clinician for a replacement.

Recommendations for all HVAD system components:

Inspect all HVAD system components for signs of damage. Damaged equipment should be reported to your clinician and replaced.

The HVAD system components were designed and tested to function for the following periods starting from when it was provided to you:

HVAD system peripherals or accessories	Expected useful life (Starting from when it was provided to you)
HVAD pump	2 years
HVAD controller (primary and back-up)	2 years
HeartWare batteries	500 charge and discharge cycles

Medtronic

Controller AC adapter	6 months
Controller DC adapter	1 year
Alarm adapter	1 year
HeartWare battery charger	1 year
Carry cases (HeartWare shoulder pack, waist pack, convertible patient pack)	1 year