Date: 30.04.2024

Urgent Field Safety Notice **Intuy Knee**

For Attention of*: Identify either by name or role who needs to be aware of the hazard and/or take action. If this is multiple recipients then include full list. Distributors of Intuy Knee

Contact details of local representative (name, e-mail, telephone, address etc.)*

This could be a distributor or local branch of the manufacturer. To be added at the appropriate stage in the different local languages

Belgium, Rebocon Bionics B.V., Sander Crevits, screvits@rbionics.com, +31640904872
Germany, Wilhelm Julius Teufel GmbH, Christoph, christoph.thiel@teufel-international.com, +49 7161 15684 137
France, Janton Solutions, Marie-Emmanuelle Ledroit, meledroit@janton-solutions.com, +33619448515
UK, Ortho Europe, Andy Hunter, andy.hunter@ortho-europe.com, +44 (0)7823 403 450
Norway, Ortopartner AS, Jorn Kvamme, jorn@ortopartner.no, +4797160767
Portugal, Padrao, Carlos Quelhas, cquelhas@padrao-ortopedico.com, + 351 22 9065149
Spain, EMO, Leopoldo Fernandez, leopoldofernandez@emo.es, +34616660665

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2 Reason for Field Safety Corrective Action (FSCA)*	
2	<p>1. Description of the product problem*</p> <p>EN: A motor encoder warning causes the knee to stop at random position.</p>

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3. Type of Action to mitigate the risk*	
3.	<p>1. Action To Be Taken by the User*</p> <p> <input type="checkbox"/> Identify Device <input type="checkbox"/> Quarantine Device <input checked="" type="checkbox"/> Return Device <input type="checkbox"/> Destroy Device </p> <p><input checked="" type="checkbox"/> On-site device modification/inspection</p> <p><input type="checkbox"/> Follow patient management recommendations</p> <p><input type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU)</p> <p><input type="checkbox"/> Other <input type="checkbox"/> None</p> <p>For "on-site device modification" the distributor must coordinate with the orthopaedic clinics to get the affected devices firmware updated.</p>
3.	<p>2. By when should the action be completed?</p> <p style="text-align: right;">Specify where critical to patient/end user safety</p> <p style="text-align: center;">31-05-2024</p>
3.	<p>3. Particular considerations for:</p> <p style="text-align: right;">Choose an item.</p> <p>Is follow-up of patients or review of patients' previous results recommended?</p> <p>Choose an item.</p> <p>Provide further details of patient-level follow-up if required or a justification why none is required</p>
3.	<p>4. Is customer Reply Required? *</p> <p><i>(If yes, form attached specifying deadline for return)</i></p> <p style="text-align: right;">No</p>
3.	<p>5. Action Being Taken by the Manufacturer</p> <p> <input type="checkbox"/> Product Removal <input type="checkbox"/> On-site device modification/inspection <input checked="" type="checkbox"/> Software upgrade <input type="checkbox"/> IFU or labelling change <input type="checkbox"/> Other <input type="checkbox"/> None </p> <p>Provide further details of the action(s) identified.</p>
3	<p>6. By when should the action be completed?</p> <p style="text-align: center;">31-05-2024</p>
3.	<p>7. Is the FSN required to be communicated to the patient /lay user?</p> <p style="text-align: right;">No</p>
3	<p>8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?</p> <p>Choose an item. Choose an item.</p>

Rev 1: April 2024

FSN Ref: Critical Software Update (Intuy Knee Firmware v2.5.10)

FSCA Ref: CAPA-38

4. General Information*	
4.	1. FSN Type* New
4.	2. For updated FSN, reference number and date of previous FSN Provide reference and date of previous FSN if relevant
4.	3. For Updated FSN, key new information as follows: Summarise any key difference in devices affected and/or action to be taken.
4.	4. Further advice or information already expected in follow-up FSN? * No
4	5. If follow-up FSN expected, what is the further advice expected to relate to: Eg patient management, device modifications etc
4	6. Anticipated timescale for follow-up FSN For provision of updated advice.
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)
	a. Company Name Rebocon Bionics B.V.
	b. Address _____
	c. Website address _____
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. * Yes
4.	9. List of attachments/appendices: EN: Appendix A. Affected Devices, Appendix B. Newsletters DE: Anhang A. Betroffene Geräte, Anhang B. Newsletter NL: Bijlage A. Betrokken apparaten, Bijlage B. Nieuwsbrieven FR: Annexe A. Appareils concernés, Annexe B. Newsletters
4.	10. Name/Signature

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback..*</p>

Rev 1: April 2024

FSN Ref: Critical Software Update (Intuy Knee Firmware v2.5.10)

FSCA Ref: CAPA-38

Note: Fields indicated by * are considered necessary for all FSNs. Others are optional.

Appendix A: Affected devices

Anhang A: Betroffene Geräte

Bijlage A: Betrokken apparaten

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EN: Below is the list of all affected devices on the European Market, the serial number can be found engraved on the side of the device.

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EN: Newsletter on 19 April:

Dear Customers,

The following software updates will be released on 19/04/2024:

- Intuy Knee Firmware **v2.5.9**

Risk

Recently we have received several reports of **W02 warning** and the knee **locks** at certain position.

There is a **risk of falling**.

The knee should not be locked and for most knees, this is a false alarm due to a **software bug**. This bug exists in Intuy Knee firmware **2.5.7 and 2.5.8**, but not **2.5.6**.

Action

We kindly urge you to update the firmware for your customers to 2.5.9 **ASAP**. After updating, please perform a power cycle.

- If the error/warnings are removed after this, the knee does not have to be returned for repair.
- If the error/warnings persist after updating the firmware and a power cycle, then please do contact our customer support to return the product for repair.

EN: Newsletter on 30 April:

Dear Customers,

The following software updates will be released on 30/04/2024:

- Intuy Knee Firmware **v2.5.10**

Risk

On 19 April, we released the firmware **v2.5.9** to prevent the knee from being **locked**.

After the update, some users experienced jerky movements from the knee. In very rare case, the knee still locks.

This is caused by the change to the software algorithm related to encoder noise and disturbance handling.

Action from Our Side

In **v2.5.10**, we **revert** the changes and made it the same as **v2.5.6**. The jerky movements are prevented after this update.

Due to the **reversion**, the following features are temporarily disabled. We shall push it out in the future release after thorough testing.

1. In 2.5.9, "Battery indication without the app" is added. In 2.5.10, it is removed.
2. In 2.5.9, we disabled "holding function" for activity level 1. In 2.5.10, holding function is enabled for activity level 1.

3. In 2.5.9, after clicking "Reset Default", Activity Level is set to 4. In 2.5.10, Reset Default sets the activity level back to 2.

Action from Your Side

We kindly urge you to update the firmware for your customers to 2.5.10 ASAP. After updating, please perform a power cycle.

- If the error/warnings are removed after this, the knee does not have to be returned for repair.
- If the error/warnings persist after updating the firmware and a power cycle, then please do contact our customer support to return the product for repair.

We are sorry for the inconvenience. We shall do our best to prevent this from happening.

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