

More information on back side.

William Cook Australia Pty. Ltd. 95 Brandl Street Eight Mile Plains, QLD 4113 Australia +61 7 38 41 11 88	William Cook Europe Aps Sandet 6 4632 Bjaeverskov Denmark +45 56 86 86 86	Cook Incorporated P.O. BOX 489 Bloomington, IN 47402 U.S.A. 800.457.4500
---	---	--



This patient has received a

Patient Name

Date of placement

Hospital

Physician name

Physician Phone#



The Cook Celect Platinum Vena Cava Filter is **MR Conditional**. A patient with this filter can be scanned safely after placement under the following conditions:

- Static magnetic field of 3.0 Tesla or 1.5 Tesla.
- Maximum spatial magnetic gradient of 1600 Gauss/cm (16.0 T/m) or less.
- Normal operating mode.
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence).

MR image quality may be compromised if the area of interest is within approximately 21 mm of the position of the Cook Celect Platinum Vena Cava Filter.

Additional information is available at
www.cookmedical.com

AL-1205-432

88,9 x 101,6 mm